

Conners' Continuous Performance Test – Second Edition (CPT-II, V.5)

This assessment is a computer-administered measure of motoric inhibitory control which, unlike conventional CPT's, requires the subject to respond to all targets except for the target letter. Thus, the motor system is primed to respond, and for a minority of the trials the subject must inhibit the prepotent response. The test is well-normed for children and adults. It provides measures of omission and commission errors, mean and SD reaction times (RT's) for hits, and a measure of between-block RT changes. It pulls strongly for errors of impulse control. It is the second edition of the instrument.

Link for purchase:

<http://www.mhs.com/product.aspx?gr=edu&prod=cptii5&id=overview>

Conners, C. K., & Staff, M. H. S. (2000). *Conners' Continuous Performance Test II (CPT II V. 5)*. Multi-Health Systems Inc., North Tonawanda, NY.